

BF046 DERMATO SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Cycloheximide	250 mg
Chlortetracycline	50 mg
Gentamicin	50 mg

Directions :

Rehydrate the contents of 1 vial with 5 ml of 50% v/v acetone. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) D.T.M. Agar Base (B160) .Mix well and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use :An antibiotic supplement, recommended for the selective isolation of pathogenic dermatophytes

Storage: Between 2 to 8°C.

Packing :5 vial