

BF051 OLEIC ALBUMIN SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Alkalinized oleic acid	0.050%
Albumin fraction V	5.0% in normal saline

Directions:

Warm up the refrigerated supplement to room temperature. Shake well and aseptically add the entire content (20 ml) of 1 vial of sterile Oleic Albumin Supplement to 180 ml of sterile, molten, cooled (45-50°C) Dubos Oleic Agar Base (B1016) / Dubos Oleic Broth Base (B1192) along with 5000 to 10000 units of penicillin. Mix well.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An enrichment supplement recommended for cultivation of *Mycobacterium tuberculosis*

Storage: Between 2 to 8°C.

Packing : 5 vial

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.