

TECHNICAL SHEET

<b>BF052</b>	<b>BASIC FUCHSIN</b>
<b>FORMULA</b>	
<b>Ingredients :</b>	
Basic Fuchsin	6.00 g
<b>Directions :</b>	
Dissolve the contents of 1 vial aseptically in 60 ml. of 95% alcohol. Add 4 ml of this (10%) solution to 1000 ml of Endo Agar Base(B515). Heat to boiling to dissolve the medium completely and sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.	
<b>Precautions :</b>	
1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials. 3. Caution: Basic Fuchsin is potential carcinogen and care should be taken to avoid inhalation of the powdered dye and contamination of the skin.	
<b>Use :</b>	
Basic Fuchsin is recommended for preparation of Endo Agar Base(B515) used for confirmation of the presumptive test for the members of the coliform group	
<b>Storage:</b>	
Store below 30°C.	
<b>Packing :</b>	
1 vial	

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.