BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF054

SODIUM BICARBONATE SOLUTION

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Sodium bicarbonate (10%w / v)

Directions:

Warm up the refrigerated solution to room temperature. Add content of 1 vial aseptically with 980 ml sterile, molten, cooled (45-50°C) GTC Agar Base (B1050) along with 2 vials of GTC Supplement (BF063). Mix well and pour into sterile petri plates.

20ml

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use :A selective solution used in media is recommended for the selective isolation of Enterococcus species.

Storage: Between 2 to 8°C.

Packing: 5 vial