BIOMARK Laboratories-INDIA www.biomarklabs.com

TECHNICAL SHEET

	CDC CUDDI EMENT
3F057	GBS SUPPLEMENT

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Metronidazole5mgGentamicin1mg

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and add it to 475 ml sterile, molten, cooled (45-50°C) GBS Medium Base (B1048) along with 25 ml inactivated Horse Serum. Mix well

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation and presumptive identification of group B Streptococci.

Storage: Between 2 to 8°C.

Packing:5 vial