BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

BF061

LINCO T SUPPLEMENT

(lincomycin-colistin-amphotericin b-trimethoprim)

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

0.500 mg
3.000 mg
0.500 mg
3.250 mg

Directions :

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix gently, avoiding frothing to dissolve and aseptically add it to 240 ml of sterile, molten, cooled (45-50°C), GC Agar Base (B1049) along with250 ml of sterile 2% w/v Haemoglobin solution (BF036). Also add sterile rehydrated contents of 1 vial of Yeast Autolysate Supplement (BF038). Mix gently and pour into sterile petri plates.

Precautions :

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement recommended for the selective isolation and cultivation of Neisseria species.

Storage: Between 2 to 8°C.

Packing :5 vial