BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF066 POTASSIUM CHLORATE SUPPLEMENT

FORMULA

Ingredients:

(Per vial)

Potassium chlorate 1.00 g Distilled water 10.00 ml

Directions:

Aseptically add the contents of 1 vial to 988 ml of sterile, molten ITC Broth Base (B1068I) along with 2 ml Ticarcillin Supplement (BF065)/2.5 ml for Yersinia ITC Broth Base (B1472) along with one vial of reconstituted Ticarcillin Irgasan Supplement (BF160) for 250ml of medium. Mix well and dispense into sterile tubes or flasks.

Precautions: 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

Use: A liquid Potassium Chlorate Supplement, recommended for the selective

enrichment of Yersinia enterocolitica.

Storage: Between 2 to 8°C.

Packing: 10 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01

Rev: December 2020