

BIOMARK Laboratories-INDIA
www.biomarklabs.com
TECHNICAL SHEET

BF067

BROMO CRESOL PURPLE

FORMULA

Ingredients:

(Per vial, sufficient for 1000 ml medium)

Bromo cresol purple 15 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile 0.2 N sodium hydroxide and add to 1000 ml of sterile, molten KF Streptococcal Agar Base (B1077)/ along with 1 vial of 1% TTC Solution (BF044). Mix Well and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: A dye supplement which can be used for identification of faecal Streptococci.

Storage: Between 10 to 30°C.

Packing : 5 vials