

BF068

KIMMIG SELECTIVE SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Part A

Colistin sulphate

40mg

Part B

Novobiocin

50mg

Directions :

Rehydrate the contents of one vial each of Part A and Part B, separately with 2.5 ml sterile distilled water and mix well. On reconstitution aseptically add to 500 ml sterile, molten, cooled (45-50°C) Kimmig Fungi Agar Base (B1084). Mix well and pour into sterile petri plates

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use :An antimicrobial supplement recommended for selective isolation and cultivation of fungi.

Storage: Between 2 to 8°C.

Packing :5 vial