

BF069

TRICHOMONAS SELECTIVE SUPPLEMENT I

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Streptomycin	500 mg
Penicillin G	125000 Unit

Directions :

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and avoid frothing of solution. Aseptically add the contents to 475 ml of sterile, molten, cooled (45-50°C) Kupferberg Trichomonas Broth Base (Trichomonas Broth Base, Kupferberg) (B1091) along with 25 ml sterile beef or human serum. Mix well and dispense as desired.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antimicrobial supplement recommended for selective isolation of Trichomonas species from clinical specimens.

Storage: Between 2 to 8°C.

Packing : 5 vial