

**BF070**

**LACTIC SUPPLEMENT**

**FORMULA**

**Ingredient :**

(Per vial sufficient for 500 ml medium)

Polysorbate 80	5 ml
Cycloheximide	3.50 mg
Phenyl ethanol	1.50 g

**Directions :**

Aseptically add the contents of 1 vial to 500 ml sterile, molten, cooled (45-50°C) Lactic Acid Bacteria Selective Agar Base B574 / Lactic Acid Bacteria Selective Broth Base (Raka ray no.3) B1251. Mix well and dispense as desired.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** A selective supplement recommended for the isolation of lactic acid bacteria in beer and brewing process.

**Storage:** Store at 2-8°C.

**Packing:** 5 Vials