

TECHNICAL SHEET

BF070

LACTIC SUPPLEMENT

FORMULA

Ingredient :

(Per vial sufficient for 500 ml medium)

Polysorbate 80	5 ml
Cycloheximide	3.50 mg
Phenyl ethanol	1.50 g

Directions :

Aseptically add the contents of 1 vial to 500 ml sterile, molten, cooled (45-50°C) Lactic Acid Bacteria Selective Agar Base B574 / Lactic Acid Bacteria Selective Broth Base (Raka ray no.3) B1251. Mix well and dispense as desired.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: A selective supplement recommended for the isolation of lactic acid bacteria in beer and brewing process.

Storage: Store at 2-8°C.

Packing: 5 Vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.