BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF073 OXFORD LISTERIA SUPPLEMENT, MODIFIED

FORMULA

Ingredients:

(per vial, sufficient for 500 ml medium)

Amphotericin B5.00 mgColistin sulphate10.00 mgAcriflavine hydrochloride2.50 mgCefotetan1.00 mgFosfomycin5.00 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix gently to dissolve and aseptically add to 500 ml Sterile, molten Listeria Oxford Medium Base (B889)/ Listeria Oxford Agar Base (B1416). Mix well and pour into sterile Petri plates.

Precautions: 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

Use: An antimicrobial supplement recommended for selective isolation of

Listeria species.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01

Rev: December 2024