

BF078

ROSOLIC ACID

FORMULA

Ingredients:

Rosolic acid 0.1 g per vial

Directions:

Dissolve the contents of one vial in 10 ml of 0.2 N NaOH Mix well and aseptically add it to 990 ml M-FC Agar (B871), M-FC Agar Base, Modified (B870), M-FC Broth (B1100). Boil to dissolve the medium completely.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: Rosolic acid solution is recommended for selective isolation of coliform bacteria.

Storage: Between 2 to 8°C.

Packing: 5 vials