

BF084

MUELLER TELLURITE SERUM

FORMULA

Ingredients :

(Per vial sufficient for 975 ml medium)

| | |
|-------------------------|------------|
| Potassium tellurite | 100.000 mg |
| Calcium D pantothenate | 0.050 mg |
| Horse serum | 12.500 ml |
| Sodium lactate solution | 10ml |
| Ethanol | 2.500 ml |

Directions:

Warm up the refrigerated solution to room temperature. Add the contents of 1 vial to 975 ml of sterile, molten, cooled (45-50°C) Mueller Tellurite Agar Base (B566). Mix well and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: This supplement is recommended for isolation, cultivation and differentiation of *Corynebacterium diphtheria*

Storage: Store below -20°C.

Packing : 5 vial

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.