

**BF085**

**MYCOPLASMA ENRICHMENT SUPPLEMENT**

**FORMULA**

**Ingredients:**

(Per vial, sufficient for 70 ml medium)

Horse serum	20 ml
Yeast extract (25%w/v solution)	10ml
Thallos acetate	25 mg
Penicillin G	20000IU

**Directions:**

Aseptically add the contents of 1 vial to 70 ml of sterile, molten, cooled (45-50°C) Mycoplasma Agar Base (B265) / Mycoplasma Broth Base w/o CV(B267)/ Mycoplasma Broth Base w/ CV (B266) PPLO Modified Broth Base w/o CVFor Mycoplasma Broth Base w/ CV (B266) add 0.285 ml of 1% Potassium tellurite (BF010) along with Mycoplasma Enrichment Supplement (BF085). Mix well and dispense as desired

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use: A** selective supplement for the isolation of Mycoplasma

**Storage:**Between 2 to 8°C.

**Packing :**5 vial

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.