BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

BF087	PARK AND SANDERS SELECTIVE SUPPLEMENT I
FORMUL	Α
Ingredie	nts:
	(Per vial, sufficient for 950 ml medium)
	Vancomycin 10.00 mg Trimethoprim 10.00 mg
Direction	Rehydrate the contents of one vial aseptically with 5 ml of 50% aqueous ethanol. Mix well and aseptically add it to 950 ml of sterile, cooled (45- 50°C) Park and Sanders Broth Base (B463)/ Park and Sanders Broth Base (B1480) along with 50 ml of sterile defibrinated lysed horse blood and Campylobacter Cefex Agar Base (B965) containing 50 ml sterile defibrinated sheep blood. Mix well and dispense.
Precautio	 I. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.
Use :	An antibiotic supplement recommended for the selective isolation and enumeration of thermo tolerant Campylobacter
Storage:	Between 2 to 8°C.
Packing :	5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01