

BF090 S.F.P. SUPPLEMENT (PERFRINGENS S.F.P. SUPPLEMENT)

FORMULA

Ingredients:

(Per vial sufficient for 500 ml medium)

Kanamycin sulphate	6mg
Polymyxin B	15000IU

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add it to 475 ml of sterile, molten, cooled (45-50°C) Perfringens Agar Base (T.S.C. / S.F.P. Agar Base, (B282)/ S.F.P. Agar Base (B374) along with 25 ml of sterile Egg Yolk Emulsion (BF003). Mix well and pour into sterile petri plates

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement, recommended for the selective isolation of Clostridium perfringens

Storage: Between 2 to 8°C.

Packing : 5 vials