

TECHNICAL SHEET

BF091 T.S.C. SUPPLEMENT (PERFRINGENS T.S.C. SUPPLEMENT)

FORMULA

Ingredients:

(per vial, sufficient for 500 ml medium)

D-Cycloserine	200.00 mg
---------------	-----------

Directions:

Rehydrated the content of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 475 ml of sterile, molten Perfringens agar Base (B282), along with 25 ml of Egg Yolk Emulsion or 500 ml of Perfringens Agar Base (B282)/Perfringes Agar Base (B1496)/TSC Agar Base (B1462)/Tryptose Cycloserine Dextrose Agar Base (B061) or 750 ml of Tryptose Cycloserine Azide Agar Base (B062). Mix well and pour into the sterile Petri plates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation of Clostridium perfringens.

Storage: Between 2 to 8°C.

Packing: 5 vials