BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF093 CFC SUPPLEMENT

FORMULA

Ingredients:

(per vial, sufficient for 500 ml medium)

Cetrimide 5.00 mg
Fucidin 5.00 mg
Cephaloridine 25.00 mg

Directions:

Rehydrate the contents of one vial aseptically with 1 ml 95% Ethanol + 1 ml 0.2N NaOH solution. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Pseudomonas Agar Base (B695). Mix well and pour into sterile Petri plates.

Precautions: 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation of

Pseudomonas species.

Storage: Between 2 to 8°C.

Packing: 1 vial