

**BF096**

**DIPHTHERIA VIRULENCE SUPPLEMENT (PART A & B)**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 500 ml medium)

Part A

Horse serum 100ml

Part B

Potassium tellurite 1ml

**Directions :**

Warm up the refrigerated contents of Part B vial and aseptically add 29ml sterile distilled water. Mix thoroughly. Aseptically add warmed up (to 50°C) contents of Part A and B vials to sterile, molten, cooled (45-50°C) Tinsdale Agar Base (B334) as required. For 10 ml of M314 : 1.0 ml of Part A and 0.3 ml of Part B, is recommended. Mix well and pour into sterile petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** A selective supplement recommended for the isolation and presumptive identification of *Corynebacterium diphtheriae*

**Storage:** Between 2 to 8°C.

**Packing :** 5 vial