BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF102 YERSINIA SELECTIVE SUPPLEMENT

FORMULA

Ingredients:

(per vial sufficient for 500 ml of medium)

Cefsulodin7.50 mgTriclosan (Irgasan)2.00 mgNovobiocin1.25 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water and 1 ml of ethanol. Mix gently to dissolve the contents completely and aseptically add to 500 ml of sterile, Molten Yersinia selective agar base (B001I). Mix well and pour into sterile Petri plates.

Precautions: 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation of

Yersinia enterocolitica.

Storage: Between 2 to 8°C.

Packing: 2 vials