

BF104

SM SELECTIVE SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Tyrothricin	20.000 mg
Cycloheximide	50.000 mg
Bacitracin	50.000 mg
Captan	10.000 mg
Vancomycin	10.000 mg
Chloromycetin	0.005 mg
Penicillin G	0.001mg

Directions :

Rehydrate the contents of one vial aseptically with 10 ml of 50% acetone. Mix well, avoiding frothing to dissolve and aseptically add it to 990 ml of sterile, molten, cooled (45-50°C) SM Selective Agar Base (B373) along with 1 vial of TTC Solution 1% (BF044). Mix well and pour Into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use :An antibiotic supplement recommended for the selective isolation and cultivation of *Pseudomonas solanacearum* .

Storage: Between 2 to 8°C.

Packing :5 vial