

BF105

S.T.A. SELECTIVE SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Streptomycin sulphate	250 mg
Cycloheximide	25 mg
Thallos acetate	25 mg

Directions :

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water. Mix well and aseptically add it to 495 ml of sterile, molten, cooled (45-50°C) S.T.A. Agar Base (B369). Mix well and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : Recommended for the selective isolation of *Brochothrix thermosphacta* from meat products.

Storage: Between 2 to 8°C.

Packing : 5 vial