## **BIOMARK Laboratories-INDIA**

## www.biomarklabs.com

### **TECHNICAL SHEET**

## BF108

## **CAT SELECTIVE SUPPLEMENT**

## **FORMULA**

# **Ingredients:**

(Per vial, sufficient for 500 ml medium)

Cefoperazone 4 mg
Teicoplanin 2 mg
Amphotericin B 5 mg

## Directions:

Rehydrate the contents of one vial aseptically with 2 ml of sterile distilled water. Mix well and add aseptically to 500 ml of sterile, molten, cooled (45-50°C) Blood Free Campylobacter Selectivity Agar Base (B945)/ Mix well and pour into sterile petri plates.

## **Precautions:**

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** An antibiotic supplement recommended for selective isolation of thermophilic Campylobacter species.

Storage: Between 2 to 8°C.

Packing: 5 vials