

BF112 CAMPYLOBACTER SELECTIVE SUPPLEMENT (KARMALI),MODIFIED

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Sodium pyruvate	50.000 mg
Vancomycin	10.000 mg
Cefoperazone	16.000 mg
Amphotericin B	5.000 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 95% ethanol. Mix well and aseptically add it to 500 ml. Sterile molten Karmali Campylobacter Agar Base (B1080) Mix well and pour into sterile Petriplates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: This supplement recommended for the selective isolation of thermo tolerant Campylobacter species.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.