

**BF114                      GEORGE KIMMIG SELECTIVE SUPPLEMENT**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 500 ml medium)

Penicillin G	20000 IU
Streptomycin	20 mg

**Directions :**

Rehydrate the content of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C), Kimmig Fungi Agar Base (B1084). Mix well and pour into sterile petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** An antimicrobial supplement recommended for selective isolation and cultivation of fungi

**Storage:** Between 2 to 8°C.

**Packing :** 5 vial