BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

BF117	FRASER SELECTIVE SUPPLEMENT
FORMULA	
Ingredients:	
(per vial sufficier	nt for 1000 ml medium)
Acriflavin hydrocl Nalidixic acid	hloride 25.00 mg 20.00 mg
Directions: Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add to 990 ml of sterile, cooled (45-50°C) Fraser Secondary Enrichment Broth Base (B1045). Mix well and dispense as desired.	
Precautions :	 For Laboratory Use. Follow proper, established laboratory procedures in handling and disposing of infectious materials.
Use :	An antimicrobial supplement recommended for selective isolation and cultivation of Listeria monocytogenes from food and environmental specimens.
Storage:	Between 2 to 8°C.
Packing :	5 vials