

BF119

M-CP SELECTIVE SUPPLEMENT I

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

D-Cycloserine	200mg
Polymyxin B sulphate	12.500mg

Directions :

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 485 ml sterile, molten, cooled (45-50°C) M-CP Agar Base (B1211) along with one vial of M-CP Selective Supplement I (BF119) / M-CP Selective Supplement II, Modified (BF120). Mix well and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement recommended for the selective isolation of *Clostridium perfringens*

Storage: Between 2 to 8°C.

Packing : 5 vial