## BIOMARK Laboratories-INDIA www.biomarklabs.com

**TECHNICAL SHEET** 

BF130 KANAMYCIN SULPHATE SELECTIVE SUPPLEMENT

**FORMULA** 

Ingredients:

(Per vial, sufficient for 500 ml medium)

Kanamycin sulphate 10mg

**Directions:** 

Add 2 ml of sterile distilled water to one vial and mix gently to dissolve completely. Add the contents to 500ml of either Kanamycin Esculin Azide Agar Base (B1133) or Kanamycin Esculin Azide Broth Base (B1145).

**Precautions:** 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

**Use:** An antimicrobial supplement recommended for the selective isolation of

group D streptococci.

**Storage:** Between 2 to 8°C.

**Packing:** 5 vial