BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF132 LEGIONELLA (GVPC) SELECTIVE SUPPLEMENT, MODIFIED

FORMULA

Ingredients:

(Per vial sufficient for 500 ml medium)

Glycine 1.500g Vancomycin hydrochloride 0.500mg Polymyxin B sulphate 40000IU Cycloheximide40mg

Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add to 440 ml sterile, molten, cooled (45-

50°C) Legionella Agar Base (B1202) along with 50 ml of Legionella Growth Supplement BCYE (BF131). Mix well and pour into sterile petri plates

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for selective isolation of Legionella species.

Storage: Between 2 to 8°C.

Packing: 5 vial

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01

Rev: December 2020