

**BF133                      LEGIONELLA BMPA SELECTIVE SUPPLEMENT**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 80 ml medium)

Cefamandole	0.400 mg
Polymyxin B	8000 IU
Anisomycin	8 mg

**Directions :**

Rehydrate the contents of 1 vial aseptically with 5 ml of distilled water. Mix well and aseptically add to 80 ml of sterile, molten, cooled (45-50°C) Legionella Agar Base (B1202) along with 15 ml of sterile Legionella Growth Supplement (BCYE) (BF131). Mix well and pour into sterile petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** An antibiotic supplement recommended for selective isolation of Legionella Species

**Storage:** Between 2 to 8°C.

**Packing :** 5 vial

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.