

**BF137 IMRV/RV SELECTIVE SUPPLEMENT**

**FORMULA**

**Ingredients:**

(Per vial, sufficient for 500 ml medium)

Novobiocin	10.00 mg
------------	----------

**Directions:**

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water and aseptically add it to 1000 ml sterile, molten, cooled (45-50°C) Semisolid IMRV Medium Base (B1257), Modified Semisolid RV Medium Base/ Rappaport Vassiliadis Semisolid Medium Base, Modified (MSRV)(B1402)/ & 500 ml of Semisolid RV Medium Base (B1259). Mix well and pour into sterile petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** An antibiotic supplement recommended for isolation of Salmonella from Foodstuffs and other materials.

**Storage:** Between 2 to 8°C.

**Packing :** 5 vials

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.