BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

BF138 SHIGELLA SELECTIVE SUPPLEMENT

FORMULA

Ingredients:

(Per vial, sufficient for 1000 ml medium)

Novobiocin 0.550 mg

Directions:

Rehydrate the contents of 1 vial with 5 ml sterile distilled water and aseptically add it to 1000 ml of sterile, molten, cooled (45-50°C) Shigella Broth Base (B1286) Mix well and dispense as desired.

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation and cultivation of Shigella species.

Storage: Between 2 to 8°C.

Packing:5 vials

Rev: December 2020