

**TECHNICAL SHEET**

**BF147 EGG YOLK EMULSION,50%**

**FORMULA**

Ingredients  
(100 ml per vial)

Egg yolk	50ml
Sterile saline	50ml

**Directions:**

Warm up the refrigerated egg yolk emulsion to room temperature. Shake well to attain uniform emulsion. (Since on refrigeration emulsion has a tendency to form layers or small lumps). Aseptically Add 80 ml emulsion in 920 ml of sterile, molten, cooled (45-50°C) Anaerobic Egg Agar Base, S.F.P. Agar Base (B374), Liver-Veal-Agar Base, Modified OR Aseptically add 50 ml emulsion in 950 ml of sterile, molten, cooled (45-50°C) into MYP Agar.

**Precautions :** 1. For Laboratory Use.  
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** Sterile stabilized emulsion of egg yolk recommended for use in various culture media in accordance with FDA BAM 1998.

**Storage:** Between 2 to 8°C.

**Packing:** 5 vials

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.