#### **BIOMARK Laboratories-INDIA**

#### www.biomarklabs.com

#### **TECHNICAL SHEET**

## **BF160**

## TICARCILLIN IRGASAN SUPPLEMENT

#### **FORMULA**

# **Ingredients:**

(Per vial, sufficient for 250 ml medium)

Ticarcillin 0.25mg Irgasan 0.25mg

## **Directions:**

Rehydrate the content of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 250 ml of sterile, cooled (45-50°C) Yersinia ITC Broth Base (B1472).

### **Precautions:**

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** For selective isolation and enumeration of Yersinia enterocolitica from clinical

specimens, water and food samples.

Storage: Between 2 to 8°C.

Packing: 5 vial

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01

Rev: December 2020