### **BIOMARK Laboratories-INDIA**

## www.biomarklabs.com

#### **TECHNICAL SHEET**

#### **BF164**

### **GENTAMYCIN SELECTIVE SUPPLEMENT**

### **FORMULA**

# **Ingredients:**

(per vial sufficient for 1000 ml medium)

Gentamycin 20.00 mg

## **Directions:**

Rehydrate the content of 1vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add it to 1000 ml sterile, molten, cooled (45-50°C) Columbia Agar BH1161.Mix well and pour into sterile Petri plates.

**Precautions:** 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

**Use:** An antibiotic supplement used for detection of Clostridium perfringens

from pharmaceutical products.

**Storage:** Between 2 to 8°C.

**Packing:** 5 vials