

TECHNICAL SHEET

BF166 ECO157: H7 Selective Supplement

FORMULA

Ingredients:

(Per vial sufficient for 1000 ml medium)

Vancomycin	8 mg
Cefixime	0.050 mg
Cefsulodin	10 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 1000 ml of sterile, cooled to 45-50°C Buffered Peptone Water (B129). Mix well and dispense into sterile tubes.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: Recommended for selective isolation of E. coli O157:H7.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.