

**BI446 Bovine Serum Albumin Powder**

**Specifications**

**Appearance** : Off white to amber

**Clarity (2% solution)** : Clear

**pH (10% solution)** : 7.0 – 8.0

**Moisture** : NMT 5%

**Microbial Load:(CFU/g)**

Total aerobic microbial count by plate method when incubated at 30-35 °c for 48 hours

**Bacterial count** : NMT 5000

Total Yeast and mold count by plate method when incubated at 20-25 °c for 120 hours

**Yeast and mold count** :NMT 250

**Test for specific organisms**

E.coli : Absent

Salmonella species : Absent

Pseudomonas aeruginosa : Absent

Staphylococcus aureus : Absent

**Precaution**

1. For Laboratory Use
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Storage** : Between 2 to 8°C.

**Packing** : 5/25/100/500g