

**TECHNICAL SHEET**

<b>BW706</b>	<b>MILK AGAR W/CETRIMIDE (TWIN PACK)</b>					
<b>Formula</b>						
<b>Ingredients :</b>		<b>gms/lit.</b>				
<b>Part A</b>						
Skim milk powder		100.00				
<b>Part B</b>						
Peptic digest of animal tissue		2.50				
Sodium chloride		1.25				
Yeast extract		0.75				
Cetrimide		0.30				
Agar		15.00				
Final pH (at 25°C) : 7.3 ± 0.2						
<b>Directions :</b>						
Suspend 19.8 grams of Part B in 250 ml of distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 20 minutes. Suspend 100 grams of Part A in 750 ml of distilled water. Sterilize by autoclaving at 15 lbs pressure (121°C) for 5 minutes. After autoclaving cool both the parts to 50°C. Aseptically add Part A solution to Part B solution, mix well and pour into sterile Petri plates.						
<b>Principles :</b>						
Peptic digest of animal tissue, yeast extract and skim milk provide nitrogen, sulphur, vitamins and other growth nutrients. Sodium chloride maintains osmotic equilibrium. Cetrimide is a quaternary ammonium compound which inhibits a wide variety of microorganisms including Pseudomonas species other than Pseudomonas aeruginosa .						
<b>QC Tests – (I) Dehydrated Medium</b>						
Colour :		Part A : White to cream Part B : Cream to yellow				
Appearance :		Part A & B : Homogeneous Free Flowing powder				
<b>(II) Rehydrated medium</b>						
pH (post autoclaving/heating) :		7.3 ± 0.2				
Colour (post autoclaving/heating) :		Cream				
Clarity (post autoclaving/heating) :		Opalescent gel				
<b>(III) Q.C. Test Microbiological</b>						
Cultural characteristics observed after 24 - 48 hours at 35-37°C.						
MICROORGANISM (ATCC )		GROWTH		PIGMENT		
Escherichia coli (25922)		Inhibited				
Pseudomonas aeruginosa (27853)		Good – Luxuriant		Blue green		
Stenotrophomonas maltophilia (13637)		Inhibited				
<b>Precautions :</b>		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
<b>Limitations :</b>		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
<b>Use :</b>		It is recommended for detection and enumeration of Pseudomonas aeruginosa in water. It is recommended by BIS.				
<b>Storage :</b>		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
<b>Packing :</b>		500 gm. bottle				
<b>Product profile:</b>		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
<b>BW706</b>		119.8 (Part A +B) g/l	4.173 L	7.3 ± 0.2	NIL	121°C / 15 minutes

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.