

BF015 CAMPYLOBACTER SUPPLEMENT-III(SKIRROW)

FORMULA

Ingredients:

(Per vial sufficient for 500 ml medium)

Vancomycin	5.000 mg
Trimethoprim	2.500 mg
Polymyxin B sulphate	1250.000 IU

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well to dissolve avoid frothing of the solution. Aseptically add the contents to 500 ml sterile molten cooled (45-50°C) Brucella Agar Base (B125) / Brucella Broth Base (B126) / Blood Agar Base no 2 (B116) / Blood Agar Base no 2 (B117) / Columbia Blood Agar Base no 2 (B994)/Columbia Blood Agar Base (w/Agar 1%) (B148). Mix gently and pour dispense into sterile Petri Plates / tubes.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: Campylobacter supplement recommended for the selective isolation of campylobacter species

Storage: Between 2 to 8°C.

Packing: 5 vials