

BF023

LEGIONELLA SELECTIVE SUPPLEMENT II, MODIFIED

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Cephalothin	2.00 mg
Colistin sulphate	8.00 mg
Vancomycin	0.25 mg
Amphotericin B	5.00 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten cooled (45-50°C) Buffered Charcoal Yeast Extract Agar Base(B127), Modified Buffered Charcoal Yeast Extract Agar Base(B627) along with rehydrated content of Legionella selective supplement IV (BF025). Mix well and pour into sterile Petri plates

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation of Legionella species from mixed culture.

Storage: Between 2 to 8°C.

Packing: 5 vials