

BF050

DOYLE'ANTIBIOTIC SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Vancomycin	7.500mg
Trimethoprim	0.500mg
Polymyxin B sulphate	10000IU
Cycloheximide	25mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 460 ml of sterile, molten, cooled (45-50°C) Doyles Enrichment Broth Base (B1013) along with 35 ml sterile lysed horse blood. Mix well and dispense as desired.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement recommended for the selective isolation of Campylobacter species.

Storage: Between 2 to 8°C.

Packing : 5 vial