

**BF053**

**FGTC ANTIBIOTIC SUPPLEMENT**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 500 ml medium)

Gentamicin                    2.500mg

**Directions :**

Rehydrate the contents of one vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 1000 ml sterile, molten, cooled (45-50°C) FGTC Agar Base (B1032) along with 2 vials of Amylose Azure (BF111) And 1 vial of Sodium Bicarbonate Solution (BF054). Mix well before pouring in to sterile Petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :**An antibiotic supplement recommended for the selective isolation of Enterococcus species

**Storage:** Between 2 to 8°C.

**Packing :**5 vial