

**BF056 FRASER LISTERIA SUPPLEMENT**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 500 ml medium)

Ferric ammonium citrate                      250 mg

**Directions :**

Rehydrate the contents of 1 vial aseptically with 1-2 ml of sterile distilled water. Mix well. For primary enrichment aseptically add 2 vials of (BF056) to 1000 ml, sterile, cooled (45-50°C) Fraser Broth Base (B963I) along with rehydrated contents of 1 vial of Fraser Selective Supplement (BF117) or add 1 vial of (BF056) to 500 ml sterile Fraser Broth Base (B963I) along with rehydrated contents of 1 vial of Fraser Selective Supplement (BF117) for secondary enrichment. Mix well before dispensing as desired.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** Recommended by ISO Committee for the selective isolation and enumeration of *Listeria monocytogenes* from food

**Storage:** Store between 10 to 30°C.

**Packing :** 5 vial

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.