

**BF063**

**GTC SUPPLEMENT**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 500 ml medium)

Gentamicin                      1.250 mg

**Directions :**

Rehydrate the content of 1vial aseptically with 2 ml sterile distilled water and aseptically add to 500 ml sterile, molten, cooled (45-50°C) GTC Agar Base (B1050) along with 10 ml of sterile 10% Sodium Bicarbonate solution (BF054) just before use. Mix well and pour into Sterile petriplates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :**                                      An antibiotic supplement recommended for selective isolation of Enterococci.

**Storage:** Between 2 to 8°C.

**Packing :** 5 vial