

BF065 TICARCILLIN SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 1000 ml. medium)

Ticarcillin 1.00 mg

Directions :

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 998 ml Of sterile, molten ITC Broth Base (B10681) containing 10 ml Potassium Chlorate Supplement (BF066). Mix well and dispense into sterile tubes.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement recommended for the selective enrichment of *Yersinia enterocolitica*.

Storage: Between 2 to 8°C.

Packing : 5 vial