

BF076

GRUFT MYCOBACTERIAL SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 400 ml medium)

Penicillin	20000 IU
Nalidixic acid	14 mg
Ribonucleic acid	20 mg

Directions :

Rehydrate the contents of 1 vial aseptically with 2.5 ml sterile distilled water. Appropriate concentrations of the supplement can be achieved by absorbing 0.05 ml of reconstituted supplement on slants prepared from 8 ml L-J Medium (slants prepared with 150 ml of sterile L.J. Medium (B230) and 250 ml of whole egg emulsion collected aseptically).

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic and enrichment supplement recommended for the selective cultivation of Mycobacteria

Storage: Between 2 to 8°C.

Packing : 5 vial