

BF077

NOVOBIOCIN SELECTIVE SUPPLEMENT

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Novobiocin 15 mg

Directions :

Rehydrate the content of 1 vial aseptically with 5 ml of sterile distilled water and aseptically add to 1000 ml of sterile, cooled Lysine Iron Cystine Broth Base (B596). Mix well and dispense as desired.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use :An antibiotic supplement, recommended for the rapid presumptive detection of Salmonella species in foods and feed materials

Storage: Between 2 to 8°C.

Packing :5 vial