

TECHNICAL SHEET

BF081 MIDDLEBROOK ADC GROWTH SUPPLEMENT

FORMULA

Ingredients:

(per vial sufficient for 450 ml medium)

| | |
|---------------------------|----------|
| Bovine albumin fraction V | 2.50 g |
| Dextrose | 1.00 g |
| Catalase | 0.0015 g |
| Distilled water | 50.00 ml |

Directions:

Warm up the refrigerated supplement to 45-50°C. Shake well to form uniform suspension. Aseptically add 1 vial to 450 ml of sterile, molten, cooled (45-50°C) Middlebrook 7H9 Broth Base (B733) . Mix well and pour into sterile screw capped tubes.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antimicrobial supplement recommended for cultivation of Mycobacteria.

Storage: Between 2 to 8°C.

Packing: 1 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.